# WIND PROFILE CODE PROCEDURE

## After wave gauges calibration process:

## **Turning on the equipment:**

- 1. Turn on the electrical boxes, Ethernet Switch, Amplifier & NI cDAQ-9184
- 2. Turn on the DC power supply positioned on the instruments carriage and set to 5[V]
- 3. Open changeVMchangeA.vi located in D:\Experiment\VIs\Wave Gauges Programs\SubVIs
- 4. Insert IP Address (192.168.1.2 for larger linear stage) and remote port (503)
- 5. run program

## These values are saved until the electrical boxes are turned off

- \* To check if the DAQ Device is connected correctly: open NI MAX ⇒ Devices and interfaces ⇒ Network Devices ⇒ cDAQ9184-1B560CB right click on this cDAQ and select the "self-test"
- If completed successfully⇒ click on cards NI 9212, NI 9401, NI 9205 to enable them.
- If not successful⇒ right click on the cDAQ and click on "reserve chassis".

### Finding and saving the safe position above the highest wave for the airflow sensors:

- 1. Turn on the wind blower at frequency 40 [Hz]. Open ctrlSpeedBlower.vi located in D:\Experiment\VIs\Blower, insert 40 [Hz] and run the program
- 2. Move the airflow sensors close to the water surface
- 3. Open WaterSensorScopeSimple.vi located in D:\Experiment\VIs\Scope\_Simple

  This VI shows drop in voltage (~1.5-4[V]) when the sensor touches water
- 4. Move the linear stage so that the exposed tips of the water sensor are in the water. The drop of voltage will be evident on the graph/

### Make sure the Pitot tube and static pressure probe are not touching the water

- 5. Move the sensors up with a step of 3-4 [mm] and run the program for 30[sec]
  - If there is a drop of voltage, repeat step e)
  - If there is no drop of voltage, double check if you are in a safe position:
  - \* run the program for 60 [sec] and if there is a drop of voltage, repeat step e)
  - \* If there is no drop of voltage run the program for 120 [sec]
  - \* If there is no drop of voltage proceed to next step, otherwise, repeat step e)
- 6. Save this point as 0 by opening LinearStageMaster.vi located in \\D:\Experiment\VIs\Linear Stage and insert IP Address (192.168.1.2) and port (503).

Write "p 0" in the command window and run program

#### Sanity check

- Write "pr p" in command window, run program.
- Check if the position equals 0 in data out
- 7. Open WindProfile.vi located in *D:\Experiment\VIs\Wind Profile Program* Insert the following inputs:
  - 1. Wind blower frequency:(e.g. 40, 32, 25, 18 [Hz])
  - 2. Seeking for zero water level:
    - \* Sampling rate (e.g. 120 [Hz])
    - \* Measurement time (e.g. 60 [s])
    - \* Zero seek step (e.g. 3 [mm])
  - 3. Airflow sensors settings:
    - \* Sampling rate (e.g. 1600 [Hz])

- \* Measurement time (e.g. 60 [s])
- 4. Linear stage settings:
  - \* remote port 503.
  - \* IP Address 192.168.1.2
- 5. Z coordinates (relative to 0[mm])
- 6. Fetch (location of airflow sensors relative to drain holes [cm])
- 7. Choose the calibration file (located in *D:\Experiment\Flume measurements\CalibrationFiles*)
- 8. Insert the Path to Data files
- 8. Run the program.

From:

https://tsailwiki.com/ - T-SAIL WiKi

Permanent link:

https://tsailwiki.com/doku.php?id=wind\_profile\_code\_procedure

Last update: 2018/05/03 06:36



https://tsailwiki.com/ Printed on 2025/08/25 00:49